Video Article

High-throughput Fluorometric Measurement of Potential Soil Extracellular Enzyme Activities

Colin W. Bell¹, Barbara E. Fricks¹, Jennifer D. Rocca¹, Jessica M. Steinweg², Shawna K. McMahon³, Matthew D. Wallenstein¹

¹Natural Resource Ecology Laboratory, Colorado State University

Correspondence to: Matthew D. Wallenstein at matthew.wallenstein@colostate.edu

URL: http://www.jove.com/video/50961

DOI: doi:10.3791/50961

Keywords: Environmental Sciences, Issue 81, Ecological and Environmental Phenomena, Environment, Biochemistry, Environmental Microbiology, Soil Microbiology, Ecology, Eukaryota, Archaea, Bacteria, Soil extracellular enzyme activities (EEAs), fluorometric enzyme assays, substrate degradation, 4-methylumbelliferone (MUB), 7-amino-4-methylcoumarin (MUC), enzyme temperature kinetics, soil

Date Published: 11/15/2013

Citation: Bell, C.W., Fricks, B.E., Rocca, J.D., Steinweg, J.M., McMahon, S.K., Wallenstein, M.D. High-throughput Fluorometric Measurement of Potential Soil Extracellular Enzyme Activities. *J. Vis. Exp.* (81), e50961, doi:10.3791/50961 (2013).

Abstract

Microbes in soils and other environments produce extracellular enzymes to depolymerize and hydrolyze organic macromolecules so that they can be assimilated for energy and nutrients. Measuring soil microbial enzyme activity is crucial in understanding soil ecosystem functional dynamics. The general concept of the fluorescence enzyme assay is that synthetic C-, N-, or P-rich substrates bound with a fluorescent dye are added to soil samples. When intact, the labeled substrates do not fluoresce. Enzyme activity is measured as the increase in fluorescence as the fluorescent dyes are cleaved from their substrates, which allows them to fluoresce. Enzyme measurements can be expressed in units of molarity or activity. To perform this assay, soil slurries are prepared by combining soil with a pH buffer. The pH buffer (typically a 50 mM sodium acetate or 50 mM Tris buffer), is chosen for the buffer's particular acid dissociation constant (pKa) to best match the soil sample pH. The soil slurries are inoculated with a nonlimiting amount of fluorescently labeled (i.e. C-, N-, or P-rich) substrate. Using soil slurries in the assay serves to minimize limitations on enzyme and substrate diffusion. Therefore, this assay controls for differences in substrate limitation, diffusion rates, and soil pH conditions; thus detecting potential enzyme activity rates as a function of the difference in enzyme concentrations (per sample).

Fluorescence enzyme assays are typically more sensitive than spectrophotometric (*i.e.* colorimetric) assays, but can suffer from interference caused by impurities and the instability of many fluorescent compounds when exposed to light; so caution is required when handling fluorescent substrates. Likewise, this method only assesses potential enzyme activities under laboratory conditions when substrates are not limiting. Caution should be used when interpreting the data representing cross-site comparisons with differing temperatures or soil types, as *in situ* soil type and temperature can influence enzyme kinetics.

Video Link

The video component of this article can be found at http://www.jove.com/video/50961/

Introduction

Extracellular enzymes (EEs) produced by soil bacteria, fungi, and archaea are involved in innumerable biogeochemical processes, and are central to the processing, stabilization, and destabilization of soil organic matter and nutrient cycling in terrestrial ecosystems ¹. By producing EEs, soil microbes decompose and transform polymeric organic matter into smaller soluble molecules, thereby liberating previously bound micro- and macronutrients, which allows plants and microbes to assimilate available nutrients from the soil. EEs have been studied for decades, primarily by measuring their activities in laboratory assays²⁻⁴, since it is very difficult to directly detect and quantify enzymes.

Extracellular enzyme activity (EEA) is most strongly controlled by the concentration of enzymes and corresponding substrates. The abundance of different C-, N-, and P-degrading enzymes in soils is controlled by numerous factors including microbial biomass, community composition, substrate availability, microclimate, and stoichiometric demands^{5,6}. However, *in situ* EEAs within the soil environment are also affected by temperature^{7,8}, the binding of enzymes to soil clays and humic properties², and diffusion constraints ⁹, which ultimately regulate the active enzyme pool, in terms of size, substrate availability, and turnover rates¹⁰⁻¹². Therefore, acknowledging *in situ* soil conditions is critical when using laboratory enzyme assays to interpret soil microbial function across different environmental sites.

Many different classes of EEA can be quantified in laboratory assays using a variety of synthetic substrates (please refer to "List of Reagents Table" for more detail). Some protocols utilize substrates in assays that are coupled to a colorimetric reaction that can be detected with a spectrophotometer, while others, including the protocol we describe here, utilize substrates that are bound to a fluorescent moiety. Fluorescence EE assays are typically more sensitive (by an order of magnitude) than colorimetric assays (which use a chromogenic moiety linked with a

²Biosciences Division, Oak Ridge National Laboratory

³Department of Bioengineering, University of Colorado

synthetic substrate) ¹²⁻¹⁴. Sensitivity in EEA detection includes two aspects: one is related to the quantity of the compound of interest detected and the other related to the lowest detectable potential enzyme activity. Methods for colorimetric P-nitrophenol (pNP)-based assays can be found in past works ^{15,16}. In brief, soil (typically sieved to <2 mm and air-dried) is incubated at optimal or field-relevant temperature and pH. The rate at which the reaction product released is determined colorimetrically with a spectrophotometer ¹⁴. The higher sensitivity of the fluorescence enzyme assays is due in part to the more sensitive detection of the fluorogenic moiety separation associated with substrate degradation rather than recording absorbance after the separation of a specific chromogenic moiety at a specific wavelength. The two most commonly used synthetic fluorescent indicators are 4-methylumbelliferone (MUB)¹⁷ and 7-amino-4-methylcoumarin (MUC)^{18,19}. MUC-linked substrates are commonly associated with N-rich synthetic substrates such as proteins and/or amino acids. Fluorescent techniques were first developed for aquatic samples^{20,21}, and their application to soils requires controls for signal quenching and interference^{22,23}. Assays can either be conducted using traditional "bench top" chemistry with large volumes, or can be employed in microplate based protocols, with increased throughput but possibly higher measurement error. While there are several widely cited protocols for fluorescent detection of EEAs in soils²⁴, many labs employ subtle variations on these protocols, often inadvertently or due to differences in laboratory equipment or reagents. Seemingly small differences in the details of protocols can strongly affect measured EEAs^{25,26} and the lack of standardized enzymes makes it challenging to calibrate assays between different laboratories. Thus, there is an important need for the dissemination of detailed protocols to encourage the standardization of EEA assays.

In our protocol, soil samples are prepared by combining soil samples with a pH buffer and homogenizing with a blender. The slurries are then inoculated with a nonlimiting amount of fluorescently labeled C-, N-, or P-rich substrate, chosen depending on the specific research question of interest. Using soil slurries in the enzyme assays serves as a control to minimize substrate diffusion limitations. The fluorescent moieties are quenched till they are cleaved from their respective substrates, and thus enzyme activity can be detected as the fluorescent dye is released from the substrate by an enzyme-catalyzed reaction. The increasing fluorescence intensity through time reflects the rate of the enzyme-catalyzed reaction.

The general concept of the fluorescence enzyme assay is that synthetic substrates bound with a fluorogenic moiety (fluorescent dye), are added to soil samples ²⁷. During enzyme-catalyzed substrate degradation, the bond breaks between the fluorescent dye and the substrate. The fluorescent dye liberated from the substrate is consequently used as an indirect assessment of enzyme activity, and can be quantified using a microplate reader to detect the fluorescence intensity of the dye. In brief, fluorescence quantification is accomplished as the liberated dye emits light of one wavelength after absorbing light of a different wavelength. The fluorescence intensity is recorded by a plate reader capable of both excitation and detection. Enzyme activity can be subsequently quantified based on the known fluorescent-dye concentrations of the substrate (*i.e.* known quantities of synthetic substrate added to soil samples) along with referencing a standard dilution curve of fluorescence intensities for the specific fluorogenic moiety of the substrate used in the assay (*i.e.* 4-methylumbelliferone (MUB) or 7-amino-4-methylcoumarin (MUC)). (Please refer to the protocol section for specific details on enzyme activity quantification).

Laboratory soil enzyme assays are useful for assessing microbial community function, but there are several technical limitations that users should recognize¹⁰. Fluorescence assays can suffer from interference caused by impurities and/or instability of many fluorescent compounds when exposed to light, so caution is required when handling fluorescent substrates²⁵. Soil particulates and/or organic material in the soil slurries can also interfere with fluorescence intensities, known as the quenching effect ²⁶. Furthermore, laboratory enzyme assays only assesses potential EEAs under laboratory conditions. *In vitro* assays measure EEAs under conditions where substrate diffusion and abundance is nonlimiting. Therefore, the data provided by these assays may not be a good proxy for EEAs under *in situ* soil conditions¹⁰. Overall, enzyme activity is very useful for relative comparisons in which soil types are similar. However, when using this method to compare activities among soils that differ in physical or chemical properties, caution should be used. This is due to the fact that differences in soil type and temperature can drastically alter the state of *in situ* enzyme kinetics. Another limitation is that relatively few substrates are commercially available (compared to the natural environment). Furthermore, the synthetic substrates used for enzyme assays are relatively simple (easily soluble) may not accurately represent the soil substrates present or available *in situ*. Another factor to consider is that using soil slurries will incorporate the activity of some stabilized enzymes (*i.e.* immobilized by organic matter or clays) that may not be active under *in situ* conditions ². Laboratory enzyme assays also do not provide information regarding the persistence of enzymes in the soil (enzyme turnover rates) or information regarding the specific microbial species that are producing soil enzymes.

Protocol

1. Assay Set-up

- 1. Label 3 deep well plates: "Sample", "MUB standard", and "MUC standard".
- 2. Pour each standard (MUB of 0, 2.5, 5, 10, 25, 50, and 100 μM) and substrate (BG, CB, NAG, PHOS, XYL, AG, LAP) into separate clean, prelabeled reservoirs (oriented in rows).
- 3. Pipette the appropriate MUB standard into corresponding wells of MUB standard plate (Table 1).
 - 1. Pipette 200 μl of 0 μM MUB into row A of MUB standard plate.
 - 2. Pipette 200 µl of 2.5 µM MUB into row B of MUB standard plate.
 - 3. Pipette 200 μ I of 5 μ M MUB into row C of MUB standard plate.
 - 4. Pipette 200 μ I of 10 μ M MUB into row D of MUB standard plate.
 - 5. Pipette 200 μ I of 25 μ M MUB into row E of MUB standard plate.
 - 6. Pipette 200 μl of 50 μM MUB into row F of MUB standard plate.
 - 7. Pipette 200 µl of 100 µM MUB into row G of MUB standard plate.
 - 8. Repeat steps 1.3.1 1.3.8 for MUC standards into corresponding wells of the MUC standard plate.

- 4. Prep soil slurries for each soil sample.
 - 1. Weigh 2.75 g of field moist soil into blender.
 - 2. Add 91 ml of 50 mM buffer. Blend on high for 1 min.
 - 3. Pour blender contents into a clean glass bowl at least as wide as 8-channel pipette with a stir bar.
 - 4. Place on stir plate, mix soil slurry on low.
 - 5. Pipette 800 µl of soil slurry into first column of the Sample plate (**Table 2**).
 - 6. Repeat into 1st column of MUB standard and MUC standard plates (**Table 1**).
 - 7. Rinse blender, stir plate and stir bar with DI H₂O or buffer between soil samples.
 - 8. Repeat steps 1.4.1 1.4.7) for each sample, moving to the next column with each subsequent soil sample (Tables 1 and 2).
- 5. Pipette the appropriate substrate (in this example, 200 µM concentrations) into corresponding wells of the Sample plate (Table 2).
 - 1. Pipette 200 µl of BG substrate into row A of Sample plate.
 - 2. Pipette 200 µl of CB substrate into row B of Sample plate.
 - 3. Pipette 200 µl of NAG substrate into row C of Sample plate.
 - 4. Pipette 200 µl of PHOS substrate into row D of Sample plate.
 - 5. Pipette 200 µl of XYL substrate into row E of Sample plate.
 - 6. Pipette 200 µl of AG substrate into row F of Sample plate.
 - 7. Pipette 200 µl of LAP substrate into row G of Sample plate.
 - 8. Repeat steps 1.5.1 1.5.7 for each incubation temperature sample plate.

2. Incubation of Assay Plates

- 1. Seal the deep-well plates ("Sample", "MUB standard", and "MUC standard") with plate mats.
- 2. Invert each (sealed) plate by hand until the solution is thoroughly mixed (~10x).
- 3. Place plates in appropriate incubator for required incubation time period (Table 3).
- 4. Record initial time as Time 0. Also record the appropriate incubation times so you know when to remove the plate from the incubator (**Table 3**).
- 5. At the end of each incubation period, centrifuge the sealed plates for 3 min at \sim 2,900 x g.
- 6. After certification, transfer 250 µl from each well of the incubated deep well plates into corresponding wells in black flat-bottomed 96-well plates (one black plate will be used for each incubated deep well plate). It is important to transfer samples from the incubated deep well plates into the same wells (i.e. A1...A12, etc.) in the black flat-bottomed 96-well plates (**Table 1**, 2).

3. Fluorescent Measurements on a Microplate Fluorometer

- 1. Following manufacturer's instructions for the fluorometric plate reader used, set the excitation wavelength to 365 nm and emission wavelength to 450 nm (or use appropriate filters).
- 2. Read the three plates on the fluorometer. Important: samples and standard plates (i.e. MUB or MUC) must be read using the same gain.
 - 1. Set the spectrophotometer to automatic gain and read the MUC and MUB Standard plates.
 - 2. Set gain to manual, and decrease the value from optimal to the next lowest number, rounded to the nearest five, for each Standard plate. Repeat 2x for each plate.
 - 3. Set the gain to automatic and run the Sample plate. Rerun the Sample plate manually to match the highest gain for each of the Standard plates (*i.e.* MUB and MUC).

4. Data Analysis

- 1. Enter standard curve fluorescence data into spreadsheet for MUB (Table 4a) and MUC standard (Table 4b)
 - 1. Convert (µM) concentrations to (µmol) (**Tables 4a** and **4b**).
- 2. For each sample's standard curve fluorescence data, calculate the slope, y-intercept, and R² values for MUB and MUC (μmol) standard concentrations. Acceptable R² values should exceed 0.98 for each sample (**Tables 4a** and **4b**). It may be useful to visualize standard curves in a scatter plot; with fluorescence reads plotted as the dependent variable (y-axis) and standard concentration (μmol) plotted as the independent variable (x-axis) (**Figure 1**).
- 3. You will use the MUB and MUC standard curve slope and y-intercept values to calculate the sample + substrate raw fluorescence data into potential EEAs. You must first enter the sample + substrate raw fluorescence data into a new spreadsheet (Table 5a). You will also need to enter the incubation time (hr) and soil dry weight for each sample (Table 5b). Note, the time-values entered in this example are 3 hr, because the samples were incubated at 25 °C (Table 3).
 - 1. Subtract the standard curve intercept values from the corresponding sample fluorescence values and then divide by the slope of the corresponding standard curve (**Table 5c**). To accomplish this, rearrange standard line equation to solve for sample enzyme concentration (x); x=(y-b)/m where: [y = sample fluorescence raw read; b = intercept from MUB or MUC Standard curve; m = slope from MUB or MUC Standard curve].
 - 2. Multiply sample µmol, from step 4.3.1 above, by 91 ml (buffer volume used in soil slurry) (Table 5d).
 - 3. Divide value obtained in step 4.3.2 by the sample-specific incubation time and dry soil mass (**Table 5e**): [enzyme amt. (µmol) x 91 ml x incubation (hr) x dry soil (g) x 0.8 ml x 1,000]

Note: 0.8 ml is the volume of slurry used in each well, and 1,000 corrects for the actual soil amount in each well.

4. Multiply the value obtained in step 4.3.3 by 1,000 to get the desired units (nmol activity/g dry soil/hr) (Table 5f).

Representative Results

Enzyme assays can be used to quantify potential EEAs and to compare activities among similar samples. Here, we present representative results from an experimental climate study comparing soils that experienced ambient climate conditions (ACN) to soils that were exposed to elevated CO₂ and heating treatments (EHN) (**Figures 2-6**). Plant cover in all plots were characteristic of a northern mixed grass prairie dominated by the C4 grass *Bouteloua gracilis* (H.B.K) Lag. and two C3 grasses, *Hesperostipa comata* Trin and Rupr. and *Pascopyrum smithii* (Rydb.); about 20% of the vegetation is composed of sedges and forbs. More information regarding the site description and field experimental design can be found in references²⁸⁻³⁰. Soils were collected from two different depths (0-5 cm and 5-15 cm) within each of the treatment plots using a 1.5 cm diameter core to assess soil EEAs in response to altered climate conditions. In our examples (*i.e.* **Figures 2-6**), the sample size is (N = 3). Regardless of this minimal sample size, the variation is relatively small in most cases (as demonstrated by the error bars) and robustly reflects variability in potential EEAs among treatment plots. Analysis of Variance (ANOVA) and Tukey post hoc multiple comparisons were used to identify significant shifts in enzyme activity among treatment plots and soil depths.

The following representative results have been provided to demonstrate how this high-throughput, fluorometric assay can be used to test (1) overall EEA in soils, (2) how EEA stoichiometry can be indicative of ecosystem-level processes and (3) the relationship between incubation temperature and EEA. Soil EEA's are commonly studied to relate shifts in microbial function to soil nutrient cycling; useful indicators to assess microbial nutrient demands in response to climate change, plant community shifts, and more broadly ecosystem functioning³¹⁻³³. EEA stoichiometry has been more recently adopted as an index to assess soil biochemical nutrient cycling by intersecting ecological stoichiometric theory and metabolic theory of ecology to assess potential microbial nutrient imbalances corresponding to environmental conditions⁵. Numerous studies have suggested that wide stoichiometric ratios are indicative of nutrient growth limitations³⁴⁻³⁶; and as soil nutrients become limited, microbes respond by allocating metabolic resources to produce specific enzymes to acquire deficient nutrients³⁷. Ecoenzymatic C:N:P stoichiometry ratios are thus useful to identify relative shifts in potential microbial community nutrient demands in response to various environmental perturbations⁵. Lastly, temperature sensitivities of EEAs can be useful to assess how soil microbial community functional diversity is likely influenced by temperature shifts^{7,38}. Enzyme temperature sensitivities can widely vary between soils for a single enzyme class, and microbial communities producing enzymes have demonstrated shifts in enzyme activity corresponding to shifts in climate from historic conditions⁷. Thus enzyme activity questions related to the thermal ecology of EEs can be a useful way to assess microbial functional dynamics and belowground ecosystem processes in response to climate changes^{39,40}.

In this example, potential C-, N-, and P-enzyme acquisition activities assayed at 0-5 cm soil depths did not differ by experimental treatment (**Figure 2a**). However, at 5-15 cm soil depths, several potential EEAs did differ significantly (**Figure 2b**). For example, the C-degrading enzymes β -1,4-glucosidase and β -D-cellobiohydrolase were lower in the EHN plots ($p \le 0.038$; **Figure 2b**) compared to the ACN plots. The N and P mineralizing enzymes (β -1,4-N-acetylglucosaminidase and phosphatase, respectively) were also lower in the EHN plots ($p \le 0.012$; **Figure 2b**) compared to ACN at 5-15 cm soil depths.

Calculating and plotting the sum of all C-, N-, or P- cycling potential EEAs can be a useful approach to observe broader patterns regarding potential soil C-, N-, and/or P-cycles (**Figure 3**). In this example, the sum of β -1,4-glucosidase, β -D-cellobiohydrolase, β -Xylosidase, and α -1,4-glucosidase potential EEAs was calculated to represent potential C cycling activities. The sum of β -1,4-N-acetylglucosaminidase and L-leucine aminopeptidase was calculated to represent potential N cycling activities. Phosphatase was used to represent potential P cycling activities. In this example, potential EEAs for total C-, N- and P-cycling trended lower in the EHN plots compared to the ACN plots at the 5-15 cm soil depths (**Figure 3**). However, this trend was only significant for total N and P cycling activities ($p \le 0.046$; **Figure 3**). Soil EEAs did not significantly differ among the treatment plots at 0-5 cm soil depths (**Figure 3**). The findings for this example suggest contrasting trends in enzyme-activity functional group (*i.e.* C-, N-, or P-degrading enzymes) among the treatment plots (ACN vs. EHN) in response to soil depth. For example, C-, N- and P-degrading soil EEAs in the ambient plots (ACN) trended higher at lower depths compared to soils exposed to elevated CO₂ and heating (EHN), which demonstrated an inverse pattern (**Figure 3a-c**).

Enzyme stoichiometry is another useful approach to assess potential enzyme activities in the environment (**Figure 4**). The microbial demand for nutrients is determined by the elemental stoichiometry of microbial biomass in relation to environmental nutrient availability ³². Likewise, microbes produce specific enzymes (*i.e.* C-, N-, or P-degrading enzymes) to meet nutrient demands within their soil environments, also referred to as ecological stoichiometry ⁴¹. The ratio of potential EEAs is one way to assess microbial nutrient demands. For example, a 1:1 ratio between two enzyme functional groups (for example in C:N nutrient acquisition) would suggest that the demand for N is high relative to the demand for C when considering microbial biomass C:N ratios at the community level is typically 8:1⁴². In this example, soil enzyme stoichiometry C:N, C:P, or N:P activities significantly differ among the treatment plots at 0-5 cm soil depths (**Figures 4a-c**). However, potential enzyme C:P and N:P ratios were higher in the EHN plots compared to the ACN plots at the 5-15 cm soil depths (p = 0.05; **Figures 4b** and **4c**). This observation suggests that there is relatively higher P mineralization EEAs compared to C and N EEA the ACN plots (compared to EHN) at the 5-15 cm soil depths. Enzyme C:N acquisition activity ratios demonstrated a similar trend to that demonstrated by total C EEAs with higher C:N due to higher potential C EEAs in the ACN plots at the lower depth compared to EHN (**Figure 4a**).

Temperature can strongly influence soil EEAs. Yet, in typical lab assays, soil enzymes are measured at a single temperature that may not correspond to *in situ* temperature conditions. Our fluorescence enzyme assay method allows us to consider *in situ* temperature effects by incorporating multiple laboratory incubation temperatures for comparison. Using laboratory incubations at multiple temperatures allows us to analyze temperature-dependent enzyme kinetics data using Arrhenius plot and Q₁₀ calculations. Arrhenius plot are used to visualize activation energy; and are plotted using the logarithm of enzyme activity (y-axis dependent variable) as a function of the inverse temperature converted to degrees Kelvin (1/K) on the x-axis (*i.e.* independent variable; **Figures 5a-c**). Activation energy is commonly defined as the minimum energy

required to catalyze a chemical reaction (*i.e.* degrade a given substrate into smaller products). For our purposes, activation energy serves as a proxy for the temperature sensitivity of enzyme catalyzed reactions. Higher activation energy indicates enzyme temperature sensitivity. Likewise, activation energies (*i.e.* **Figures 5d-f**) directly correspond to Q_{10} values (*i.e.* **Figures 6a-c**). A further explanation of formulas used to calculate activation energy and practical application can be found in many past works^{40,43-45}. Arrhenius plots, activation energy, and Q_{10} plots provide redundant information and should not all be used in the same manuscript to present data for publication (**Figures 5** and **6**). Therefore, when using these techniques, it is necessary to choose the most appropriate plot type for your enzyme temperature - kinetics data¹⁰. All plot types (Arrhenius and Q_{10}) were presented here for demonstration purposes; to provide visual examples of how to present enzyme results.

In our examples, we assessed potential enzyme kinetics for potential C-, N-, and P-EEAs among both treatment plots at the two soil depths (**Figure 5**; **Figure 6**). The finding demonstrated that the temperature sensitivity of EEAs was not significantly different among treatment plots at either soil depths for activation energy as demonstrated in the Arrhenius plots (**Figures 5a-c**) and corresponding activation energies (**Figure d-f**) and (not surprisingly) for Q₁₀ plots (in this example between the 15 °C and 25 °C laboratory incubations; **Figures 6a-c**). This suggests that enzyme kinetics were not influenced by the field experimental treatments in this particular study.

	1	2	3	4	5	6	7	8	9	10	11	12
Г	S1 + MUB 0µM	S2+ MUB 0µM	S3+ MUB 0µM	S4+ MUB 0µM	S5+ MUB 0µM	S6+ MUB 0µM	S7+ MUB 0µM	S8+ MUB 0µM	S9+ MUB 0µM	S10+ MUB 0µM	S11+MUB 0µM	S12+MUB 0µM
	\$1+ MUB 2.5µM	S2+MUB 2.5μM	53+MUB 2.5µM	S4+MUB 2.5µM	S5+MUB 2.5μM	56+MUB 2.5µM	\$7+MUB 2.5µM	S8+MUB 2.5µM	S9+MUB 2.5µM	\$10+MUB 2.5µM	511+MUB 2.5µM	\$12+MUB 2.5µM
	S1+ MUB 5µM	S2+ MUB 5µM	S3+MUB SµM	S4+ MUB 5µM	S5+ MUB 5µM	S6+ MUB 5µM	S7+ MUB SµM	S8+ MUB 5µM	S9+ MUB 5µM	S10+ MUB 5µM	S11+MUB SµM	S12+MUB 5µM
	S1+MUB 10µM	S2+MUB 10µM	S3+MUB 10µM	S4+MUB 10µM	S5+MUB 10µM	S6+MUB 10µM	57+MUB 10µM	S8+MUB 10µM	S9+MUB 10µM	\$10+MUB 10µM	S11+MUB 10µM	\$12+MUB 10µM
Γ	S1+ MUB 25µM	S2+ MUB 25μM	\$3+ MUB 25µM	S4+ MUB 25µM	S5+ MUB 25µM	S6+ MUB 25µM	S7+ MUB 25µM	S8+ MUB 25µM	59+ MUB 25µM	\$10+ MUB 25µM	511+ MUB 25µM	\$12+ MUB 25µM
Г	S1+ MUB 50µM	S2+ MUB 50μM	\$3+ MUB 50µM	S4+ MUB 50µM	S5+ MUB 50µM	S6+ MUB 50µM	57+ MUB 50µM	S8+ MUB 50µM	S9+ MUB 50µM	510+ MUB 50µM	511+ MUB 50µM	S12+ MUB 50µM
	S1+ MUB 100µM	S2+ MUB 100μM	31+MUB 100μM	S4+ MUB 100μM	S5+ MUB 100μM	56+ MUB 100μM	87+ MUB 100µM	S8+ MUB 100μM	59+ MUB 100µM	\$10+MUB 100µM	511+MUB 100μM	512+MUB 100µM
Г												

Table 1. Deep well plate design for fluorescence standard concentrations with soil samples. Template for organizing and loading fluorescence standards (MUB or MUC) with soil samples into the deep well plate prior to incubating. Note: The columns represent the same soil samples (800 μl of soil slurry additions). A gradient of fluorescence standard concentrations is loaded for each row (200 μl additions). Each cell represents fluorescence standard concentration plus soil slurry additions. For example, S1 - S12 represent soil samples (800 μl of soil slurry); MUB 0-100 μM = 4-Methylumbelliferone concentrations (200 μl additions). In this example, row H is open, and is consequently available to include a higher fluorescence standard concentration if needed. This decision to increase the standard curve concentrations may be necessary for higher potential enzyme activity (and/or the substrate concentrations used) for your specific soil samples. The potential enzyme activity for your samples should fall within the range of the standard curve values. Click here to view larger table.

1	2	3	4	5	6	7	8	9	10	11	12
S1 + (BG)	S2 + (BG)	S3 + (BG)	S4 + (BG)	S5 + (BG)	S6 + (BG)	S7 + (BG)	S8 + (BG)	S9 + (BG)	S10 + (BG)	S11 + (BG)	S12 + (BG)
S1 + (CB)	1.87 ± 0.00	S3 + (CB)	S4 + (CB)	S5 + (CB)	S6 + (CB)	S7 + (CB)	S8 + (CB)	S9 + (CB)	S10 + (CB)	S11 + (CB)	S12 + (CB)
SI +	100000000000000000000000000000000000000	S3 +	S4 +	S5+	S6 +	S7 +	S8+	S9 +	S10 +	S11 +	S12 +
(NAG		(NAG)	(NAG)	(NAG)	(NAG)	(NAG)	(NAG)	(NAG)	(NAG)	(NAG)	(NAG)
S1 +	and the second second	S3 +	S4 +	S5 +	S6 +	S7 +	S8 +	S9 +	S10+	S11 +	S12 +
(PHOS		(PHOS)	(PHOS)	(PHOS	(PHOS)	(PHOS)	(PHOS)	(PHOS)	(PHOS)	(PHOS)	(PHOS)
S1 +	The second second	S3+	S4+	S5+	S6+	S7+	S8+	S9+	S10+	S11 +	S12+
(XYL		(XYL)	(XYL)	(XYL)	(XYL)	(XYL)	(XYL)	(XYL)	(XYL)	(XYL)	(XYL)
S1 + (A	G) S2 + (AG)	S3 + (AG)	S4 + (AG)	S5 + (AG)	S6 + (AG)	S7 + (AG)	S8 + (AG)	S9 + (AG)	S10 + (AG)	S11 + (AG)	S12+ (AG)
S1 +	S2 +	S3 +	S4+	S5+	S6+	S7+	S8+	S9+	S10 +	S11 +	S12 +
(LAP	(LAP)	(LAP)	(LAP)	(LAP)	(LAP)	(LAP)	(LAP)	(LAP)	(LAP)	(LAP)	(LAP)

Table 2. Deep well plate design for soil samples with substrate. Template for organizing and loading soil samples and substrates into the deep well plate prior to incubating. Note: The columns represent the same soil samples (800 μl of soil slurry additions). Different substrates are loaded across each row (200 μl additions). Each cell represents soil samples plus substrate addition. For example, S1 - S12 represent unique soil samples (800 μl of soil slurry). Substrates (200 μl additions) include: BG = 4-Methylumbelliferyl β-D-glucopyranoside; CB = 4-Methylumbelliferyl β-D-cellobioside; NAG = 4-Methylumbelliferyl N-acetyl-β-D-glucosaminide; PHOS = 4-Methylumbelliferyl phosphate: XYL = 4-Methylumbelliferyl-β-D-xylopyranoside; AG = 4-Methylumbelliferyl α-D-glucopyranoside; and LAP = L-Leucine-7-amido-4-methylcoumarin hydrochloride. In this example, row H is open, and is consequently available to include another substrate to assess another potential enzyme activity if desired. Click here to view larger table.

Temperature	Time
4 °C	~23 hr
15 °C	6 hr
25 °C	3 hr
35 °C	1.5 hr

Table 3. Incubator temperatures required for corresponding incubation time periods. Incubation temperatures and corresponding time periods for the enzyme assay procedure.

		sample id	51	52	53	54	55	S6	87	58	59	510	S11	512
		Treatment	ACN	ACN	ACN	ACN	ACN	ACN	DID	EHN	HIEN	EBIN	EHON	DID
	meta data:	DIP	3	4	5	3	4	5	3	4	5	3	4	5
		depth	0.5	0.5	0.5	5 to 15	5 to 15	5 to 15	0.5	0.5	0.5	5 to 15	5 to 15	5 to 15
		well #	1	2	3	4	5	6	7	В	9	10	11	12
MUB (uM)	MUB (umols)							MUB umels	17.00	376.3	100		100	70
0	0		63	48	252	34	41	48	122	55	65	71	62	34
2.5	0.0005	190-0110	349	320	1278	349	358	428	491	962	291	1330	1451	301
5	0.000	MUB (raw	674	587	2258	637	662	782	831	1706	519	2524	2740	554
10	0.002	fluorescence	1280	1082	4415	1285	1341	1495	1576	3398	1019	4943	5434	1056
25	0.005	data):	3214	2802	11901	3525	3535	4162	4315	7846	2505	13455	13578	2966
50	0.01		6302	6018	22158	6807	6697	7932	8087	17263	5392	24135	24632	5752
100	0.02		10564	10250	36893	12121	12366	15285	14122	29399	9092	39383	39298	10246
		slope:	534773.94	521232.34	1869135.01	613273.16	621626.37	766373.63	711527.47	1492507.06	461949.76	2000736.26	1993503.92	518326.53
		intercept:	265.31	148.51	1027.61	163.85	152.48	89.52	307.17	452.50	156.85	1258.95	1492.16	136,20
		R-squared:	0.99	0.993	0.989	0.996	0.998	0.999	0.994	0.993	0.993	0.985	0.982	0.996
MUC [µM]	MUC (µmols)) MUC umols						
0	0		29	29	79	30	54	39	58	28	23	64	47	27
2.5	0.0005		669	562	1509	1262	1156	1130	1142	655	517	1209	1242	1344
5	0.000	MUC (now	1345	1378	2501	2297	2100	2131	2474	1339	1151	2339	2346	1508
10	0.002	fluorescence	2483	2115	4712	4723	4389	4415	4635	2657	1920	4647	4411	5090
25	0.005	data):	6429	5282	12395	11654	11058	10832	11787	6401	4972	11205	11039	12710
50	0.01	100000	12322	10590	23607	22821	21705	21493	23095	12773	9176	21989	21174	24935
100	0.02	5 666	23493	19735	42659	41071	39402	41674	43103	25582	18754	40263	39131	44162
		slope:	1176249.61	987750.39	2143998.43	2069246.47	1986135.01	2086872.84	2165284.14	1276178.96	929758.24	2021897.96	1961378.34	2244056.5
		intercept:	212.06	237.52	702.58	598.86	485.40	195.63	418.65	43.16	102.47	553.28	553.85	482.83
		R-squared:	0.999	0.999	0.997	0.997	0.997	1.000	0.999	1.000	1.000	0.998	0.998	0.995

Table 4. MUB and MUC standard curve calculations. Demonstrates how to organize a) MUB and b) MUC raw fluorescence data (μmol) to subsequently calculate slope, y-intercept, and R-squared values. To calculate slope, y-intercept, and R-squared values from your standard concentration curves, select (or plot) the MUB and/or MUC raw fluorescence data as the dependent variable (y-axis) and standard concentration (μmol) as the independent variable (x-axis). Note: MUB = 4-Methylumbelliferone; MUC = 7-Amino-4-methylcoumarin. Click here to view larger

	sumple id	51	52	53	54	35	56	57	55	59	510	511	512	
	Treat	ACN	ACN	ACN	ACN	ACN	ACN	EHIN	EHN	EHN	EHN	EHIN	EHN	
etadata:	rep	3	4	5	3	4	5	3	4	5	3	4	5	
	dopth	0.5	0.5	0.5	5 to 15	5 to 15	5 10 15	0.5	0.5	0.5	5 10 15	5 10 15	5 10 15	
	well #	1	2	3	4	5	6	7	8	9	10	11	12	
	a) Enter raw sample	florescence data (l	belowit	0080	11053155	975/5-21	5,525	94.00	55, 7575	98552	V2010	0.085		
	80	1065	1216	2143	2514	2646	4775	1515	2521	2045	3390	3239	561	
	CB	1322	1481	2656	3291	3350	6017	1975	2009	2190	3613	2198	922	
il duny+	NAG	2264	3130	3526	4207	6410	8397	3116	3516	3905	5221	3055	1471	
ific substrate	PHOS	10306	11694	15844	20403	20074	24323	9560	17419	20568	29265	16490	4643	
fluorescence	XYL	1358	1626	2199	3201	2445	8526	3306	1765	2425	3992	2636	1065	
data)	AG	702	863	1855	2064	1628	5443	1763	1417	1080	2781	1755	487	
	LAP	15996	15779	23440	27237	20579	22917	15785	16915	14598	27768	33112	12737	
	b) Exter data (helow): incubation time(h) and soil dry wt													
	incubation (tracib)	3	3	3	3	3	3	3	3	3	3	3	3	
	soil dry wt	2.402	2.263	2.462	2.547	2,456	2,654	2.558	2.604	2,416	2.809	2.604	2.848	
	S) Calculation: raw data values subtracted intercept and divided by the depe													
	BG	0.001495	0.002048	0.000597	0.004321	0.004011	0.006114	0.002123	0.001386	0.004094	0.00106	0.000872	0.00083	
	CB	0.001976	0.002556	0.000671	0.005099	0.003144	0.007734	0.002344	0.001043	0.004401	0.001177	0.000354	0.00151	
	NAG	0.003737	0.00572	0.001337	0.006593	0.010066	0.01084	0.003948	0.002254	0.007897	0.00198	0.000784	0.00257	
	PHOS	0.018776	0.02215	0.007927	0.033902	0.032047	0.001621	0.013004	0.011368	0.044183	0.013995	0.007923	0.00569	
	XYL	0.002043	0.002835	0.000625	0.004952	0.003688	0.011008	0.004215	0.000879	0.00491	0.001366	0.000574	0.00179	
	AG	0.000817	0.001371	0.000443	0.003095	0.002374	0.000983	0.002046	0.000848	0.001995	0.000781	0.000132	0.00067	
	LAP	0.029416	0.029988	0.011991	0.044145	0.042512	0.029786	0.021753	0.01103	0.031261	0.01325	0.015861	0.02431	
	d) Calculation: provious value multiplied by 91 (b): the sturry is made with 91mL)													
	BG	0.136079	0.186369	0.054303	0.399239	0.365026	0.556359	0.199226	0.126119	0.37254	0.096472	0.07933	0.07457	
	CB	0.179811	0.232635	0.079279	0.464019	0.468065	0.703836	0.213309	0.094902	0.400913	0.10707	0.002221	0.13795	
	NAG	0.340107	0.520528	0.121635	0.599939	0.916039	0.986439	0.359232	0.205077	0.718654	0.180207	0.071341	0.23434	
	PHOS	1.708577	2.015684	0.721345	3,003168	2.916314	2,877509	1.18338	1.034489	4.020616	1.273906	0.684626	0.79123	
	XYL	0.185937	0.25795	0.056835	0.450664	0.335602	1.001757	0.383532	0.080025	0.446806	0.124308	0.052214	0.16306	
	AG	0.074309	0.124741	0.040252	0.281991	0.218001	0.635678	0.156192	0.025907	0.151523	0.009228	0.011998	0.06155	
	LAP	2.676818	2.728869	1.091161	4.017225	3,868583	2.710559	1.979519	1.003739	2.844779	1.205718	1.443391	2.21229	
	6 Calculation: previous administration time * odd dry weight * (Le. O.Beal. used in each well) is much units: (used activity) dry solits).													
	BG	0.023605	0.094315	0.00919	0.06433	0.061928	0.087346	0.031474	0.02018	0.064249	0.01431	0.012694	0.01091	
	CB	0.001191	0.042833	0.013417	0.075909	0.079412	0.110499	0.034745	0.013183	0.069073	0.015882	0.003135	0.02015	
	NAG	0.058997	0.09584	0.020585	0.098145	0.155408	0.154867	0.058514	0.092814	0.12394	0.026731	0.011415	0.03428	
	PHOS	0.296381	0.371131	0.12208	0.491292	0.49476	0.451757	0.192798	0.165526	0.893435	0.188947	0.109547	0.11575	
	XYL	0.032254	0.047494	0.009619	0.073725	0.056996	0.157272	0.062473	0.012805	0.077057	0.018439	0.008355	0.02385	
	AG	0.01289	0.022967	0.006817	0.048125	0.035845	0.099799	0.030328	0.00941	0.031363	0.010289	0.00192	0.0090	
	LAP	0.464338	0.502443	0.184667	0.657182	0.656315	0.425546	0.322439	0.160609	0.490614	0.178847	0.230957	0.32365	
	f) Calculation; multi-	ply by 1000 to not	amal activity's d	lry soli'h	10000									
	BG	23,60511	343146	9,19023	64,3305	61,92765	87,34601	31,47411	20,18036	64,24879	14,30999	12.69354	10.9110	
	CB	31.19123	42,83304	13,4171	75,90935	79.41182	110.4992	34.74479	15.18528	69.07301	15.88196	5.155612	20.1835	
	NAG	58.99719	95,84024	20.58547	98.14464	155,4084	154,8667	58.51448	32.81441	123.94	26,73058	11.4153	34,2848	
	PHOS	298.3309	371,1305	122,0798	491,2917	494,7802	451,7566	192,758	165,5256	693,4354	188,9472	109.5471	115,759	
	XYL	32.25388	47,49407	9,618685	73,72466	56,93572	157,2715	62,47262	12.80481	77,05675	18,43894	8.354847	23,8565	
	AG	12.89007	22.96739	0.817251	40.12475	30.64514	99,79978	30.32833	9.409714	31.36238	10.26874	1.919855	9.00034	
	LAP	464.3383	502,4431	184,6671	657,1825	656,3149	425.5462	322,4393	160,6085	490.6144	178,8474	230,9573	323,656	
	Total C	99,9403	147.6091	39.04326	200,0883	234,9203	454,9133	159,0098	57,59016	241.7411	58,89964	28.12385	63,9614	
		523,3355	598,2834	205,2526	755,3271	811,7233	580.4129	380.9538	199,423	614,5544	205,578	242,3727	357,941	
	Total N													

Table 5. Enzyme activity calculations. Demonstrates how to a) organize sample raw fluorescence data and subsequently perform the steps needed (b-f) to calculate EEAs into: μmol activity/g dry soil/hr. Note: S1 - S12 represent unique soil samples. Substrates include: BG = 4-Methylumbelliferyl β-D-glucopyranoside; CB = 4-Methylumbelliferyl β-D-cellobioside; NAG = 4-Methylumbelliferyl N-acetyl-β-D-glucosaminide; PHOS = 4-Methylumbelliferyl phosphate: XYL = 4-Methylumbelliferyl-β-D-xylopyranoside; AG = 4-Methylumbelliferyl α-D-glucopyranoside; and LAP = L-Leucine-7-amido-4-methylcoumarin hydrochloride. Click here to view larger table.




Figure 1. MUB standard curve plot. Scatterplot visualization of standard curves with raw fluorescence data as the dependent variable (y-axis) and standard concentration (µmol) as the independent variable (x-axis).




Figure 2. C, N and P cycling enzyme activities. Potential EEAs at the Prairie Heating and Elevated CO_2 Enrichment (PHACE) site in control plots (ACN: exposed to ambient climate conditions) and treatment plots (EHN: exposed to increased temperature and elevated CO_2). Soil enzymes were assayed at a) 0-5 cm soil depths and b) 5-15 cm soil depths. Vertical bars represent mean \pm S.E. Note: ACN = ambient - climate plots; EHN = elevated CO_2 and heating plots. Click here to view larger figure.




Figure 3. Total C, N and P cycling enzyme stoichiometric ratios. The sum of potential EEAs involved in the acquisition of a) C, b) N, and c) P- for both treatment groups in the PHACE experiment at 0-5 cm and 5-15 cm soil depths. Vertical bars represent mean \pm S.E. Note: ACN = ambient - climate plots; EHN = elevated CO₂ and heating plots. Click here to view larger figure.




Figure 4. Enzyme stoichiometry for total C, N, and P cycling enzyme activities. Enzyme acquisition activity stoichiometric ratios at the Prairie Heating and Elevated CO_2 Enrichment (PHACE) site in control plots (ACN: exposed to ambient climate conditions) and treatment plots (EHN: exposed to increased temperature and elevated CO_2). Total soil A0 C:N, A0 C:P and A0 N:P was plotted for both treatment groups at 0-5 cm and 5-15 cm soil depths. Vertical bars represent mean A1 S.E. Note: ACN = ambient - climate plots; EHN = elevated A2 and heating plots. Click here to view larger figure.




Figure 5. Total C, N and P cycling enzyme activation energies. Temperature sensitivity of potential EEAs displayed using Arrhenius plots for total $\bf a$) C, $\bf b$) N, and $\bf c$) P degrading enzymes; as well as the corresponding activation energy values for total $\bf d$) C, $\bf e$) N, and $\bf f$) P degrading enzymes among treatment plots and soil depth at the Prairie Heating and Elevated CO₂ Enrichment (PHACE) site. Vertical bars for activation energy (*i.e.* $\bf d$ - $\bf f$) represent mean $\bf t$ S.E. Note: ACN = ambient - climate plots; EHN = elevated CO₂ and heating plots. For publication, it is important to note that the author would typically choose to present either Arrhenius plots (*i.e.* $\bf a$ - $\bf c$) or activation energy values (*i.e.* $\bf d$ - $\bf f$), but not both plot- types. Click here to view larger figure.




Figure 6. Total C, N and P cycling enzyme Q10 (15-25 °C). Temperature sensitivities of potential EEAs measured as Q_{10} , based on laboratory incubations at 15 °C and 25 °C for total a) C, b) N, and c) P degrading enzymes among treatment plots and soil depth at the Prairie Heating and Elevated CO_2 Enrichment (PHACE) site. Vertical bars for Q10 represent mean \pm S.E. Note: ACN = ambient - climate plots; EHN = elevated CO_2 and heating plots. Click here to view larger figure.

Discussion

The laboratory-based measurement of potential soil EEAs can provide important insights into microbial responses to their abiotic environment, and their consequences for ecosystem functioning. The results of this example data set suggest that minimal differences exist in soil enzyme activity or kinetics among climate treatment plots. However, inverse trends among plots encourage further investigation of covariates that may influence the production of microbial EEAs such as soil moisture, soil pH or plant growth. Overall, assessing EEAs in terms of (1) overall EEA in soils, (2) EEA stoichiometry (3) Arrhenius plots/activation energy, and (4) Q_{10} provides a wide spectrum of approaches that can be indicative of ecosystem-level processes from which to robustly characterize soil ecosystem functional dynamics.

High-throughput fluorescence-based EEA assays are a useful tool that is widely used to examine potential EEAs in soils and other environments. Importantly, potential activities reflect the enzyme pool size, but do not by themselves quantify enzyme production or turnover rates⁴⁶. Although the technique is relatively straightforward, seemingly minor differences among lab protocols can hamper the comparability of results ¹³. Unfortunately, we do not currently have suitable standardized positive controls for EEAs. The use of stoichiometric ratios is one approach for overcoming these challenges. Otherwise, the advent of high-throughput techniques has advanced the study of enzymes in the environment⁴. Careful interpretation of the data generated by these assays may elucidate important trends in microbial activity.

The robustness of the protocol stems from the ability to select for conditions specific to individual samples, but this can also result in limitations. A series of modifications will be required to ensure samples are accurately measured for individual field sites:

Buffer

The buffer you select will depend on the pH of the soil. Buffering also stabilizes the fluorescence intensity of the fluorescent standards, which is highly pH dependent^{27,47}. The pH buffer that we typically use to make the soil slurry is a 50 mM sodium acetate or Tris buffer which is chosen for the buffer's particular acid dissociation constant (pKa) to best match soil - sample pH levels. Sodium acetate has a pKa of 4.76, and Tris has a pKa of 8.06 so the amounts of these two buffers will vary in order to reach the desired pKa for the individual sample. Phosphate buffer (pKa = 7.2) has been suggested for neutral/slightly basic soils. However, we caution to test for analytical variability in preliminary studies before using this buffer, as high phosphate concentrations may interfere with enzyme activity.

Handling and storage of fluorescent substrates

Fluorescence assays can suffer from interference caused by impurities and/or instability of many fluorescent compounds when exposed to light, so caution is required when handling fluorescent substrates. We strongly recommend minimizing any light exposure to fluorescently labeled substrates and MUB and MUC standards. Using amber glass bottles or covering the glassware and containers used for making and storing fluorescent substrates and standards is highly recommended; aluminum foil to wrap glassware and containers works well. Likewise, efficiently inoculating plates and transferring to dark incubators is best practice. We recommend storing substrates and standards (-20 °C) for no longer than two months (while protecting them from light); and thawing substrates (5 °C) ~24-48 hr prior to beginning the enzyme assay(s).

Design and Replication

To best account for well to well sample variations, implementing negative assay controls and (if possible) assay replicates is recommended. Variation typically occurs due to differences in the amount of soil particles in each well and pipetting errors. Therefore, strong mixing and good pipetting technique will substantially minimize well to well variation. Furthermore, we strongly suggest implementing a negative assay control (buffer + substrate solution) to monitor substrate inconsistencies over time. This can be easily monitored when reading the enzyme plates by comparing negative control wells (we typically use the last column on the 96-well plates). Negative substrate controls are typically stable, therefore the signal increases well above detection limit, it is indicative of contamination or substrate instability, requiring replacement of the substrate and/or standard solutions.

To maximize throughput, our protocol includes several potential EEAs in a single deep well microplate, although other protocols perform one type of assay per plate (*i.e.* one substrate per plate) since different EEAs occur at different rates. Regardless, enzyme optimization must be performed on soils before performing this high throughout approach or the single plate approach. Multiple substrates can be used on a single plate if the reaction rates are relatively consistent for each enzyme within the time period of the incubation.

Our protocol recommends ~2.75 g soil to make the slurry solution, while ~1 gram is recommended in other protocols ²⁴. We suggest that using more soil (if possible) is an effective approach to better capture within-sample soil variation in enzyme activity levels. In this protocol, we incubate 800 µl of soil slurry with 200 µl substrate, while others only incubate 200 µl of soil slurry with 50 µl of substrates. This is simply a function of scaling up that ultimately does not change the measured activity. There are also practical advantages for using larger volumes. One, it is easier to avoid soil particulates while pipetting into corresponding black flat-bottomed 96-well plates before recording intensities on the fluorometer. Second, the additional volume is useful in the case of accidental spills while transferring the black flat-bottomed 96-well plates to the fluorometer before recording enzyme-related fluorescence intensities. Even slight deviations in volume will significantly decrease fluorescence among wells. Lastly, due to the high throughput nature of this protocol, we commonly elect to rely on experimental replication to represent variation in enzyme activity levels rather than performing assay replicates. It is always best practice to include analytical replicates, but in practice, we feel our protocol provides a balanced tradeoff between analytical and experimental replicates given limited resources. Likewise, our approach uses relatively more well homogenized soil (2.75 g) per assay²⁵, which inherently decreases within-soil variations. The decision to use analytical replicates should be carefully considered by testing for analytical error in preliminary studies⁴⁸. However, we recommend that experimental designs with fewer than 4 treatment-group replicates should strongly consider using assay replicates.

Soil, Buffer Volumes, and Substrate concentration optimization

The soil buffer:substrate concentration ratio is an important variable which strongly influences the measured fluorescence when performing enzyme assays. The amount of soil in the slurry added to assay or the concentration of substrate may need to be adjusted depending on the activity of enzymes in the soil sample. The amounts selected for this example were based on previous testing of these soils to ensure that substrate availability was nonlimiting, and that we were measuring maximum potential rates under the assay conditions (V_{max}) ^{44,45}. However, for soil samples that have high concentrations of enzymes, the amount of soil or substrate concentration must be increased. We have found that increasing substrate concentrations (and adjusting the standard curve accordingly) can affect linearity of the curve. Therefore, we recommend reducing soil amounts and/or proportional adjustments to buffer volumes as necessary. Regardless, it is important to optimize substrate concentration for each soil type assayed because measured EEAs can differ by an order of magnitude greater between saturated and subsaturated substrate concentrations ²⁶. Therefore differences among experimental treatments (*etc.*) are susceptible to type II statistical error and less likely to be detected at sub-saturating conditions, and statistical error ^{27,35}. To optimize substrate concentrations, preliminary enzyme assays should performed on representative soil samples using a wide range of substrate concentrations. After recording the fluorescence, simply plot the data (similarly to the standard curve example; **Figure 1**) to identify the substrate concentration (x-axis) that corresponds to enzyme activity (y-axis) where the slope levels off (~0). Likewise, the substrate concentration corresponding to the point where the slope levels off (~0) is a good indication of the optimal substrate concentration for that particular soil.

This assay ultimately measures fluorescence over a given time produced by the fluorogenic moiety that is cleaved from the substrate as a result of enzyme-mediated substrate depolymerization. Therefore, step 5 (substrate addition) is critical and must be performed as efficiently as possible

in order to minimize the time between when the substrate is added to the soils and when the assays are incubated. Likewise, once the substrate comes in contact with the sample, enzymatic reactions will start to occur. We recommend using a multi-channel pipette for this reason. It is strongly suggested to become efficient using the multi-channel pipette prior to the day you are performing your enzyme assays. To accomplish this, you can practice pipetting with water until you can easily transfer volumes into the 96-well plates with your pipettor.

Soil Quenching

Quenching refers to decreasing fluorescence intensity caused by particulates and/or organic material in the soil slurries-incubations²⁶. Quenching can be influenced by adjusting the soil:buffer ratios²⁵. Due to background fluorescence from individual samples, it is critical to run standards with samples to account for background (quenching) fluorescence. Although some protocols use a single concentration after testing that the signal is linear, we strongly recommend implementing a standard quench control for each sample to best control for the effects of quenching. Failure to do so will result in a standard curve not applicable to the sample, and an incorrect estimation of enzymatic activity. The addition of standards to soil slurries is not time-sensitive, as the standard addition does not affect the background fluorescence of the sample.

NaOH addition

The addition of NaOH is used in some protocols to optimize fluorometric enzyme activity measurements because the fluorescent dye released from the synthetic substrates exhibits peak fluorescence at pH>9.0^{26,49}. When considering this suggestion, the NaOH concentration necessary to the soil slurry pH (*i.e.* to pH>9) will vary depending on the particular soil and the pH buffer used²⁶. However, others argue that NaOH may not be necessary because signal intensity is generally very high even at lower pH, and because it introduces an additional source of measurement error. For example, the effect of NaOH additions on slurry pH and thus MUB or MUC fluorescence changes over time²⁵. MUB linked substrates have been shown to demonstrate consistent increased fluorescence for 20 min following NaOH additions before levels taper; while MUC has demonstrated steady decreased fluorescence for up to 60 min²⁶. Therefore, it is important to standardize the time between NaOH addition and fluorescence measurement. Alternatively, if fluorescence levels are sufficiently detectable without addition, conducting the assays without adding NaOH has been suggested as an equally acceptable alternative²⁶.

Temperature

Temperature sensitivity should be taken into account when deciding incubation temperature. If the primary interest is understanding enzyme kinetics, using three or more temperatures as illustrated using Arrhenius plots (in the results section) is a robust approach. If the sample site has characteristically low temperature such as permafrost soil, then the duration of incubation may need to be extended to allow for the enzymes to react in the colder incubation temperatures. While traditional enzyme kinetics suggests that an increase in temperature should result in increased enzyme activity, we have found that enzymes may be site specific in terms of temperature sensitivity⁵⁰. Therefore to understand site-specific enzyme activity potential it is critical that incubation temperature and duration be adjusted to reflect field site values.

Conclusion

EEs are critical drivers of biogeochemical processes in soils, and thus we need to be able to measure their activities. There are many challenges to measuring EEAs in soils, including interference and inhibition. Despite these challenges, standardized protocols (as the one described here) can be universally applied to measure EEAs for a wide range of enzymes. While it is fairly easy to generate quality data following these protocols, the interpretation of this data in an ecological context requires careful consideration of what these assays are really measuring, and how EEAs under assay conditions may differ from those under *in situ* conditions.

Disclosures

No conflicts of interest declared.


Acknowledgements

This publication was funded by the Enzymes in the Environment Research Coordination Network Research supported by the US National Science Foundation (DEB# 1021559). This research was supported by the US National Science Foundation (DEB# 1021559), and the US Department of Energy's Office of Science (Biological and Environmental Research). Any opinions, findings, and conclusions or recommendations expressed in this material are those of the authors and do not necessarily reflect the views of the US NSF.

References

- 1. Burns, R.G. & Dick, R.P. Enzymes in the environment: activity, ecology, and applications. Vol. 86 CRC (2002).
- Burns, R.G. Enzyme activity in soil: Location and a possible role in microbial ecology. Soil Biology and Biochemistry 14, 423-427, doi:http://dx.doi.org/10.1016/0038-0717(82)90099-2 (1982).
- 3. Bandick, A.K. & Dick, R.P. Field management effects on soil enzyme activities. Soil Biology and Biochemistry 31, 1471-1479 (1999).
- Burns, R.G., et al. Soil enzymes in a changing environment: Current knowledge and future directions. Soil Biology and Biochemistry 58, 216-234, doi:http://dx.doi.org/10.1016/j.soilbio.2012.11.009 (2013).
- Sinsabaugh, R.L., Hill, B.H., Follstad Shah & J., J. Ecoenzymatic stoichiometry of microbial organic nutrient acquisition in soil and sediment. Nature 462, 795-798, doi:http://www.nature.com/nature/journal/v462/n7274/suppinfo/nature08632 S1.html (2009).
- Sinsabaugh, R.L., et al. Stoichiometry of soil enzyme activity at global scale. Ecology Letters 11, 1252-1264, doi:10.1111/j.1461-0248.2008.01245.x (2008).

- Wallenstein, M., Allison, S.D., Ernakovich, J., Steinweg, J.M. & Sinsabaugh, R. Vol. 22 Soil Biology eds Shukla, G. & Varma, A. 245-258 Springer; Berlin Heidelberg (2011).
- 8. Wallenstein, M.D., McMahon, S.K. & Schimel, J.P. Seasonal variation in enzyme activities and temperature sensitivities in Arctic tundra soils. *Global Change Biology* **15**, 1631-1639, doi:10.1111/j.1365-2486.2008.01819.x (2009).
- 9. Steinweg, J.M., Dukes, J.S. & Wallenstein, M.D. Modeling the effects of temperature and moisture on soil enzyme activity: Linking laboratory assays to continuous field data. *Soil Biology and Biochemistry* **55**, 85-92, doi:http://dx.doi.org/10.1016/j.soilbio.2012.06.015 (2012).
- 10. Wallenstein, M.D. & Weintraub, M.N. Emerging tools for measuring and modeling the *in situ* activity of soil extracellular enzymes. *Soil Biology* and *Biochemistry* **40**, 2098-2106, doi:10.1016/j.soilbio.2008.01.024 (2008).
- 11. Sinsabaugh, R.L. Phenol oxidase, peroxidase and organic matter dynamics of soil. *Soil Biology and Biochemistry* **42**, 391-404, doi:10.1016/j.soilbio.2009.10.014 (2010).
- 12. Allison, S. Soil minerals and humic acids alter enzyme stability: implications for ecosystem processes. *Biogeochemistry* **81**, 361-373, doi:10.1007/s10533-006-9046-2 (2006).
- 13. Deng, S., Popova, I. E., Dick, L. & Dick, R. Bench scale and microplate format assay of soil enzyme activities using spectroscopic and fluorometric approaches. *Applied Soil Ecology* **64**, 84-90, doi:http://dx.doi.org/10.1016/j.apsoil.2012.11.002 (2013).
- Sinsabaugh, R., Klug, M.J., Collins, H.P., Yeager, P.E. & Peterson, S.O. in Standard Soil Methods For Long Term Ecological Research eds Robertson, G.P., Coleman, D.C., Bledsoe, C.S., & Sollins, P. Oxford University Press (1999).
- 15. Tabatabai, M.A. in *Methods of Soil Analysis* Vol. 2 eds Page, A.L., Miller, R.H., & Keeney, D.R. 903-947 American Society of Agronomy, Inc. (1994).
- 16. Parham, J.A. & Deng, S.P. Detection, quantification and characterization of B-glucosaminidase activity in soil. *Soil Biology and Biochemistry* **32**, 1183-1190 (2000).
- 17. Jacks, T.J. & Kircher, H.W. Fluorometric assay for the hydrolytic activity of lipase using fatty acyl esters of 4-methylumbelliferone. *Analytical Biochemistry* **21**, 279-285, doi:http://dx.doi.org/10.1016/0003-2697(67)90190-X (1967).
- 18. Smith, R.E., Bissell, E.R., Mitchell, A.R. & Pearson, K.W. Direct photometric or fluorometric assay of proteinases using substrates containing 7-amino-4-trifluoromethylcoumarin. *Thrombosis Research* 17, 393-402, doi:http://dx.doi.org/10.1016/0049-3848(80)90074-2 (1980).
- Jasinski, J.P. & Woudenberg, R.C. 7-Amino-4-methylcoumarin. Acta Crystallographica Section C 50, 1954-1956, doi:doi:10.1107/ S0108270194005834 (1994).
- 20. Hoppe, H.G. SIGNIFICANCE OF EXOENZYMATIC ACTIVITIES IN THE ECOLOGY OF BRACKISH WATER MEASUREMENTS BY MEANS OF METHYLUMBELLIFERYL-SUBSTRATES. *Mar. Ecol.-Prog. Ser.* 11, 299-308, doi:10.3354/meps011299 (1983).
- 21. Somville, M. Measurement and Study of Substrate Specificity of Exoglucosidase Activity in Eutrophic Water. *Applied and Environmental Microbiology* **48**, 1181-1185 (1984).
- 22. Darrah, P.R. & Harris, P.J. A fluorimetric method for measuring the activity of soil enzymes. *Plant and Soil* **92**, 81-88, doi:10.1007/bf02372269 (1986).
- 23. Freeman, C., Liska, G., Ostle, N.J., Jones, S.E. & Lock, M.A. The use of fluorogenic substrates for measuring enzyme activity in peatlands. *Plant and Soil* 175, 147-152, doi:10.1007/bf02413020 (1995).
- 24. Saiya-Cork, K.R., Sinsabaugh, R.L. & Zak, D.R. The effects of long term nitrogen deposition on extracellular enzyme activity in an Acer saccharum forest soil. *Soil Biology and Biochemistry* **34**, 1309-1315, doi:10.1016/s0038-0717(02)00074-3 (2002).
- 25. DeForest, J.L. The influence of time, storage temperature, and substrate age on potential soil enzyme activity in acidic forest soils using MUB-linked substrates and I-DOPA. *Soil Biology and Biochemistry* **41**, 1180-1186, doi:10.1016/j.soilbio.2009.02.029 (2009).
- 26. German, D.P., et al. Optimization of hydrolytic and oxidative enzyme methods for ecosystem studies. Soil Biology and Biochemistry 43, 1387-1397, doi:10.1016/j.soilbio.2011.03.017 (2011).
- 27. Marx, M.C., Wood, M. & Jarvis, S.C. A microplate fluorimetric assay for the study of enzyme diversity in soils. *Soil Biology and Biochemistry* 33, 1633-1640, doi:http://dx.doi.org/10.1016/S0038-0717(01)00079-7 (2001).
- 28. Parton, W.J., Morgan, J.A., Wang, G. & Del Grosso, S. Projected ecosystem impact of the Prairie Heating and CO2 Enrichment experiment. *New Phytologist* **174**, 823-834, doi:10.1111/j.1469-8137.2007.02052.x (2007).
- 29. Morgan, J., et al. C4 grasses prosper as carbon dioxide eliminates desiccation in warmed semi-arid grassland. Nature 476 (2011).
- 30. Carrillo, Y., Pendall, E., Dijkstra, F.A., Morgan, J.A. & Newcomb, J., M. Response of soil organic matter pools to elevated CO₂ and warming in a semi-arid grassland. *Plant and Soil* (2011).
- 31. Allison, S.D. & Vitousek, P.M. Responses of extracellular enzymes to simple and complex nutrient inputs. *Soil Biology and Biochemistry* 37, 937-944, doi:10.1016/j.soilbio.2004.09.014 (2005).
- 32. Moorhead, D.L. & Sinsabaugh, R.L. A theoretical model of litter decay and microbial interaction. Ecological Monographs 76, 151 (2006).
- 33. Allison, S.D., Weintraub, M.N., Gartner, T.B. & Waldrop, M.P. Vol. 22 Soil Biology eds Girish Shukla & Ajit Varma) 229-243 Springer Berlin Heidelberg (2011).
- 34. Treseder, K.K. & Vitousek, P.M. Effects of soil nutrient availability on investment in acquisition of N and P in Hawaiin rain forests. *Ecology* 82, 946-954, doi:10.1890/0012-9658(2001)082[0946:eosnao]2.0.co;2 (2001).
- 35. Fujita, Y., de Ruiter, P., Wassen, M. & Heil, G. Time-dependent, species-specific effects of N:P stoichiometry on grassland plant growth. *Plant and Soil* 334, 99-112, doi:10.1007/s11104-010-0495-y (2010).
- 36. Elser, J.J., Dobberfuhl, D.R., MacKay, N.A. & Schampel, J.H. Organism Size, Life History, and N:P Stoichiometry. *BioScience* 46, 674-684 (1996).
- 37. Allison, V.J., Condron, L.M., Peltzer, D.A., Richardson, S.J. & Turner, B.L. Changes in enzyme activities and soil microbial community composition along carbon and nutrient gradients at the Franz Josef chronosequence, New Zealand. *Soil Biology and Biochemistry* 39, 1770-1781, doi:http://dx.doi.org/10.1016/j.soilbio.2007.02.006 (2007).
- 38. German, D.P., Marcelo, K.R.B., Stone, M.M. & Allison, S.D. The Michaelis Menten kinetics of soil extracellular enzymes in response to temperature: a cross-latitudinal study. *Global Change Biology* **18**, 1468-1479, doi:10.1111/j.1365-2486.2011.02615.x (2012).
- 39. German, D.P., Chacon, S.S. & Allison, S.D. Substrate concentration and enzyme allocation can affect rates of microbial decomposition. *Ecology* **92**, 1471-1480, doi:10.1890/10-2028.1 (2011).
- 40. Fierer, N., Craine, J.M., McLauchlan, K. & Schimel, J.P. Liter quality and the temperature sensitivity of decomposition *Ecology* **86**, 320-326, doi:10.1890/04-1254 (2005).
- 41. Sterner, R.W. & Elser, J.J. in *Ecological Stoichiometry: The Biology of Elements from Molecules to the Biosphere* eds R.W. Sterner & J.J. Elser) Ch. 5, 179 230 Princeton University Press; Princeton, NJ, (2002).



- 42. Cleveland, C. & Liptzin, D. C:N:P stoichiometry in soil: is there a "Redfield ratio" for the microbial biomass? *Biogeochemistry* **85**, 235-252, doi:10.1007/s10533-007-9132-0 (2007).
- 43. Logan, S.R. The origin and status of the Arrhenius equation. Journal of Chemical Education 59, 279, doi:10.1021/ed059p279 (1982).
- 44. Trasar-Cepeda, C., Gil-Sotres, F. & Leirós, M.C. Thermodynamic parameters of enzymes in grassland soils from Galicia, NW Spain. *Soil Biology and Biochemistry* **39**, 311-319, doi:10.1016/j.soilbio.2006.08.002 (2007).
- 45. Menten, L. & Michaelis, M. Die kinetik der invertinwirkung. Biochem Z 49, 333-369 (1913).
- 46. Wallenstein, M.D., Haddix, M.L., Lee, D.D., Conant, R.T. & Paul, E.A. A litter-slurry technique elucidates the key role of enzyme production and microbial dynamics in temperature sensitivity of organic matter decomposition. *Soil Biology and Biochemistry* **47**, 18-26, doi:10.1016/j.soilbio.2011.12.009 (2012).
- 47. Chrost, J. Fluorescence correction for measurements of enzyme activity in natural waters using methylumbelliferyl-substrates. *Archiv für Hydrobiologie* **106**, 79 (1986).
- 48. Reed, G.F., Lynn, F. & Meade, B.D. Use of Coefficient of Variation in Assessing Variability of Quantitative Assays. *Clinical and Diagnostic Laboratory Immunology* **9**, 1235-1239, doi:10.1128/cdli.9.6.1235-1239.2002 (2002).
- 49. Mead, J.A.R., Smith, J.N. & Williams, R.T. Studies in detoxification. 67. Biosynthesis of the glucuronides of umbelliferone and 4-Methylumbelliferone and their use in fluorimetric determination of b-glucuronidase. *Biochemical Journal* 61, 569-574 (1955).
- 50. Haddix, M.L., et al. The Role of Soil Characteristics on Temperature Sensitivity of Soil Organic Matter. Soil Sci. Soc. Am. J. 75, 56-68, doi:10.2136/sssaj2010.0118 (2011).